


Product Name: Recombinant Human ODAD4 Protein
Catalog #: PEH9372



Summary

Name	TTC25_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96NG3
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HNNEAQQAIL SALDDANKGI IRELKRTNYV ENLKEKSEGE ASLYEDRIIT REKDMRRVRD EPEKVVQWD HSEDEKETDE DDEAFGEALQ SPASGKQSVE AGKARSDLGA VAKGLSGELG TRSGETGRKL LEAGRRESRE IYRRPSGELE QRLSGEFSRQ EPEELKKLSE VGRREPEELG KTQFGEIGET

image



Background

Alternative Names Tetratricopeptide repeat protein 25 (TPR repeat protein 25)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ODAD4 Protein
Catalog #: PEH9372

