


Product Name: Recombinant Human SPAG8 Protein
Catalog #: PEH9378



Summary

Name	SPAG8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99932
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSRSRLDIQ PSSEGLGPTS EPFPSSDDSP RSALAAATAA AAAAAASAAAA TAAFTTAKAA ALSTKTPAPC SEFMEPSSDP SLLGEPCAGP GFTHNIAHGS LGFEPVYVSC IAQDTCTTDD HSSNPGPVPG

image



Background

Alternative Names	Sperm-associated antigen 8 (HSD-1) (Sperm membrane protein 1) (SMP-1) (Sperm membrane protein BS-84)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use