


Product Name: Recombinant Human AIFM2 Protein
Catalog #: PEH9382



Summary

Name	FSP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9BRQ8
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VQLLSERVS NLEELPLNEY REYIKVQTDK GTEVATNLVI LCTGIKINSS AYRKAFESRL ASSGALRVNE HLQVEGHSNV YAIGDCADVR TPKMAYLAGL HANIAVANIV NSVKQRPLQA YKPGALTFL

image



Background

Alternative Names Ferroptosis suppressor protein 1 (FSP1) (EC 1.6.5.-) (Apoptosis-inducing factor homologous mitochondrion-associated inducer of death) (AMID) (p53-responsive gene 3 protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AIFM2 Protein
Catalog #: PEH9382

