


Product Name: Recombinant Human UQCC2 Protein
Catalog #: PEH9383



Summary

Name	UQCC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BRT2
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEWVPVDETKR GRDLGAYLRQ RVAQAFREGE NTQVAEPEAC DQMYESLARL HSNYYKHKYP RPRDTSFSGL SLEEYKLILS TDTLEELKEI DKGMWKKLQE KFAPKGPEED HKA

image



Background

Alternative Names	Ubiquinol-cytochrome-c reductase complex assembly factor 2 (Breast cancer-associated protein SGA-81M) (Mitochondrial nucleoid factor 1) (Mitochondrial protein M19)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human UQCC2 Protein
Catalog #: PEH9383

