


Product Name: Recombinant Human NAA11 Protein
Catalog #: PEH9385



Summary

Name	NAA11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BSU3
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LMNMQHCNLL CLPENYQMKY YLYHGLSWPQ LSYIAEDEDG KIVGYVLAKM EEEPDDVPHG HITSLAVKRS HRLGLAQKL MDQASRAMIE NFNKYVSLH VRKSNRPALH LYSNTLNFQI SEVEPKYYAD GEDAYAMKRD LSQMADELRR QMDLKKGGYV VLGSRENQET QGSTLSDSEE ACQQKNPATE ESGSDSKEPK ESVESTNVQD

image



Background

Alternative Names	N-alpha-acetyltransferase 11 (EC 2.3.1.255) (N-terminal acetyltransferase complex ARD1 subunit homolog B) (hARD2) (NatA catalytic subunit Naa11)
--------------------------	--

Note

Product Name: Recombinant Human NAA11 Protein
Catalog #: PEH9385



For Research Use Only , Not for Diagnostic Use