


Product Name: Recombinant Human DESI2 Protein
Catalog #: PEH9386



Summary

Name	DESI2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BSY9
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VYDMYWMNEY TSSIGIGVFH SGIEVYGREF AYGGHPYPFS GIFEISPGNA SELGETFKFK EAVVLGSTDF LEDDIEKIVE ELGKEYKGNA YHLMHKNCNH FSSALSEILC GKEIPRWINR LAYFSSCIPF LQSCLPKEWL TPAALQSSVS QELQDELEEA EDAAASASVA

image



Background

Alternative Names	Deubiquitinase DESI2 (EC 3.4.19.12) (Desumoylating isopeptidase 2) (DeSI-2) (PPPDE peptidase domain-containing protein 1) (Protein FAM152A)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DESI2 Protein
Catalog #: PEH9386

