


Product Name: Recombinant Human KATNAL1 Protein
Catalog #: PEH9391



Summary

Name	KATL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BW62
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GARGPVGRAH PISKSEKPST SRDKDYRARG RDDKGRKNMQ DGASDGEMPK FDGAGYDKDL VEALERDIVS RNPSIHWDDI

image



Background

Alternative Names	Katanin p60 ATPase-containing subunit A-like 1 (Katanin p60 subunit A-like 1) (EC 5.6.1.1) (p60 katanin-like 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use