


Summary

Name	NADAP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BWU0
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LFILQGPEED REAESELTVT QLKELRKQQQ ILLEKKMLGE DSDEEEEMDT SERKINAGSQ DDEMGCTWGM GEDAVEDDAE ENPIVLEFQQ EREAFYIKDP KKALQGFFDR EGEELEYEFD EQGHSTWLCR VRLPVDDSTG KQLVAEAIHS GKKKEAMIQC SLEACRILDY LGLLRQEAVS RKRKAKNWED EDFYDSDDDT FLDRTGLIEK KRLNRMKKAG KIDEKPETF SLVAKLNDAE RELSEISERL KASSQVLSES PSQDSLDAFM SEMKSGSTLD GVSRRKHLR TFE LRKEQQR LKGLIKIVKP AEIPELKKTE TQTTGAENKA KKLTLPFGA MKGGSKFKLK TGTVGKLPK RPELPPTLMR

image



Background

Alternative Names Kanadaplin (Human lung cancer oncogene 3 protein) (HLC-3) (Kidney anion

Product Name: Recombinant Human SLC4A1AP Protein
Catalog #: PEH9392



exchanger adapter protein) (Solute carrier family 4 anion exchanger member
1 adapter protein)

Note

For Research Use Only , Not for Diagnostic Use