



Product Name: Recombinant Human RAB11FIP5 Protein
Catalog #: PEH9393



Summary

Name	RFIP5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXF6
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RKSSLTQSNT SLGSDSTLSS ASGSLAYQGP GAELLTRSPS RSSWLSTEGG RDSAQSPKLF THKRTYSDEA NQMRVAPPRA LLDLQGHLDA ASRSSLCVNG SHIYNEEPQG PVRHRSSISG SLPSSGSLQA VSSRFSEEGP RSTDDTWPRG SRSNNSSEAV LGQEELSAQA KVLAPGASHP GEEEGARLPE GKPVQVATPI VASSEAVA EK

image



Background

Alternative Names	Rab11 family-interacting protein 5 (Rab11-FIP5) (Gamma-SNAP-associated factor 1) (Gaf-1) (Phosphoprotein pp75) (Rab11-interacting protein Rip11)
--------------------------	--

Note



Product Name: Recombinant Human RAB11FIP5 Protein
Catalog #: PEH9393

For Research Use Only , Not for Diagnostic Use