


Product Name: Recombinant Human SKIC8 protein
Catalog #: PEH9401



Summary

| | |
|------------------------------|--|
| Name | WDR61_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9GZS3 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 21 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | QEQAHDDAIW SVAWGTTNKE NSETVVTGSL DDLVKVWKWR DERLDLQWSL EGHQLGVVSV DISHTLPIAA SSSLDAHIRL WDLENGKQIK SIDAGPVD TLAFSPDSQY LATGTHVGKV NIFGVESGKK EYSLDTRGKF ILSIAYSPDG KYLASGAIDG IINIFDIATG |

image



Background

| | |
|--------------------------|---|
| Alternative Names | WD repeat-containing protein 61 (Meiotic recombination REC14 protein homolog) (SKI8 homolog) (Ski8) [Cleaved into: WD repeat-containing protein 61, N-terminally processed] |
|--------------------------|---|

Note

Product Name: Recombinant Human SKIC8 protein
Catalog #: PEH9401



For Research Use Only , Not for Diagnostic Use