


Product Name: Recombinant Human UBA5 Protein
Catalog #: PEH9404



Summary

Name	UBA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9GZZ9
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QERSLQVPRS GDGGGGRVRI EKMSSEVVDS NPYSRLMALK RMGIVSDYEK IRTFAVAIVG VGGVGSVTAE MLTRCGIGKL LLFDYDKVEL ANMNRLFFQP HQAGLSKVQA AEHTLRNINP DVLFEVHNYN ITTVENFQHF MDRISNGGLE EGKPVDLVLS CVDNFEARMT INTACNELGQ TWMESGVSEN AVSGHIQLII PGESACFACA PPLVVAANID EKTLKREGVC AASLPTTMGV VAGILVQNVL KFLNFGTVS FYLGYNAMQD FFPTMSMKPN PQCDNRNCRK QQEYKKKVA ALPKQEVQIE EEEIHEDNE WGIELVSEVS EEELKNFSGP VPDLP EGITV AYTIPKKQED SVTELTVEDS

image



Background

Alternative Names	Ubiquitin-like modifier-activating enzyme 5 (Ubiquitin-activating enzyme 5)
--------------------------	---

Product Name: Recombinant Human UBA5 Protein
Catalog #: PEH9404



(ThiFP1) (UFM1-activating enzyme) (Ubiquitin-activating enzyme E1 domain-containing protein 1)

Note

For Research Use Only , Not for Diagnostic Use