


Product Name: Recombinant Human NBPF3 Protein
Catalog #: PEH9407



Summary

Name	NBPF3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H094
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MNILEINKKS RPQLAENKQQ FRNLKQKCLV TQVAYFLANR QNNYDYEDCK DLIKSMLRDE RLLTEEKLAELGQAEELRQ YKVLVHSQER ELTQLREKLQ EGRDASRSLN QHLQALLTPD EPDNSQGRDL REQLAEGCRL AQHLVQKLSP ENDDDEDEDV KVEEAEKVQE LYAPREVQKA EEKEVPEDSL EECITCSNS HHPCESNQPY GNTRITFEED QVDSTLIDSS SHDEWLDAVC IIPENESDHE QEEKGVPSP RNLQESEEEE APQESWDEGD WTLSIPPDMS ASYQSDRSTF HSVEEQVGL ALDIGRHWCD QVKKEDQEAT

image



Background

Alternative Names	Neuroblastoma breakpoint family member 3 (Protein AE2) (Protein SHIIIa4)
--------------------------	--

Product Name: Recombinant Human NBPF3 Protein
Catalog #: PEH9407



Note

For Research Use Only , Not for Diagnostic Use