


Product Name: Recombinant Human GBP3 Protein
Catalog #: PEH9409



Summary

Name	GBP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H0R5
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VQLPAETLQE LLDLHRVSER EATEVYMKNS FKDVDHLFQK KLAAQLDKKR DDFCKQNQEA SSDRCSALLQ VIFSPLEEEV KAGIYSKPGG YCLFIQKLQD LEKKYYEPR KGIQAEILQ TYLKSKEVT DAILQTDQIL TEKEKEIEVE CVKAESAQAS AKMVEEMQIK YQQMMEEKEK SYQEHVKQLT EKMERERAQL

image



Background

Alternative Names Guanylate-binding protein 3 (EC 3.6.5.-) (GTP-binding protein 3) (GBP-3) (Guanine nucleotide-binding protein 3)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GBP3 Protein
Catalog #: PEH9409

