


Product Name: Recombinant Human ISCU Protein
Catalog #: PEH9410



Summary

Name	ISCU_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H1K1
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DHYENPRNVG SLDKTSKNVG TGLVGAPACG DVMKLQIQVD EKGKIVDARF KTFGCGSAIA SSSLATEWVK GKTVEEALTI KNTDIAKELC LPPVKLHCSM LAEDAIIKAAL ADYKLLKQEPK

image



Background

Alternative Names Iron-sulfur cluster assembly enzyme ISCU, mitochondrial (NifU-like N-terminal domain-containing protein) (NifU-like protein)

Note

For Research Use Only , Not for Diagnostic Use