


Product Name: Recombinant Human MED28 Protein
Catalog #: PEH9413



Summary

Name	MED28_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H204
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSTLVDELES SFEACFASLV SQDYVNGTDQ EEIRTGVDQC IQKFLDIARQ TECFFLQKRL QLSVQKPEQV IKEDVSELRN ELQRKDALVQ

image



Background

Alternative Names	Mediator of RNA polymerase II transcription subunit 28 (Endothelial-derived protein 1) (Mediator complex subunit 28) (Merlin and Grb2-interacting cytoskeletal protein) (Magain) (Tumor angiogenesis marker EG-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use