


Product Name: Recombinant Human CCDC68 Protein
Catalog #: PEH9414



Summary

Name	CCD68_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2F9
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AQRLFENYQT QSEEVRRKKQE DSKQLLQV NK LEKEQKLKQH VENLNQVAEK LEEKHSQITE LENLVQRMEK EKRTLLEKRL SLENKLLQLK

image



Background

Alternative Names Coiled-coil domain-containing protein 68 (Cutaneous T-cell lymphoma-associated antigen se57-1) (CTCL-associated antigen se57-1)

Note

For Research Use Only , Not for Diagnostic Use