


Product Name: Recombinant Human ACTL8 Protein
Catalog #: PEH9418



Summary

Name	ACTL8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H568
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TETPLREPAD RKKMLEILFE LLHVPSVLLA DQLQMSLYAS GLLTGVVVDS GYGLTRVQPF HQGRPLPASG KTLEFAGQDL SAYLLKSLFK EDCRRCLFQ LETVAVTQMN KCYVPQNLGE ALDFRERQQS ALDESNTYQL PDGSRVELTP MQRVAPEMFF SPQVFEQPGP SIPRAIVESV ESCEISLRPL LVSHVMACGG

image



Background

Alternative Names Actin-like protein 8 (Cancer/testis antigen 57) (CT57)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ACTL8 Protein
Catalog #: PEH9418

