


Product Name: Recombinant Human HAUS4 Protein
Catalog #: PEH9420



Summary

Name	HAUS4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H6D7
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTNVTSEDKK FHETLEQRLL VTELMRLLGP SQEREIPPLL GLEKADLLEL MPLSEDFVWM RARLQQEEVEE QLKKKCFTLL CYYDPNSDAD SETVKAARKVW KLAEVLVGEQ QQCQDAKSQQ

image



Background

Alternative Names HAUS augmin-like complex subunit 4

Note

For Research Use Only , Not for Diagnostic Use