


Product Name: Recombinant Human EXO5 Protein
Catalog #: PEH9423



Summary

| | |
|------------------------------|--|
| Name | EXO5_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9H790 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 42 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | SAEASGFSDL SDSEFLEFLD LEDAQESKAL VNMPGPSSSES LGKDDKPISL QNWKRGLDIL SPMERFHLKY LYVTDLATQN WCELQTAYGK ELPGFLAPEK AAVLDTGASI HLARELELHD LVTVPVTTKE DAWAIKFLNI LLLIPTLQSE GHIREFPVFG EGEVLLVGV IDELHYTAKG ELELAELKTR RRPMLPLEAQ KKKDCFQVSL YKYIFDAMVQ GKVTPASLIH HTKLCLEKPL GPSVLRHAQQ GGFSVKSLGD LMELVFLSLT LSDLPVIDIL KIEYIHQETA TVLGTEIVAF KEKEVRAKVQ HYMAYWMGHR EPQGV DVEEA WKCRTCTYAD ICEWRKGGV LSSTLAPQVK |

image



Background

Alternative Names Exonuclease V (Exo V) (hExo5) (EC 3.1.-.-) (Defects in morphology protein 1)

Product Name: Recombinant Human EXO5 Protein
Catalog #: PEH9423



homolog)

Note

For Research Use Only , Not for Diagnostic Use