


Product Name: Recombinant Human PHAX Protein
Catalog #: PEH9424



Summary

Name	PHAX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H814
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TVAPSDRPLQ LPKVLGGDSA MRAFQNTATA CAPVSHYRAV ESVDSEESF SDSDDDSCLW KRKRQKCFNP PPKPEPFQFG QSSQKPPVAG GKKINNIWGA VLQEQNQDAV ATELGILGME GTIDRSRQSE TYNILLAKKL RKESQEHTKD LDKELDEYMH GGKKMGSKKEE ENGQGHLKRK RPVKDRLGNR PEMNYKGRIE ITAEDSQEKV ADEISFRLQE PPKDLIARVV RIIGNKKAIE LLMETAEEVEQ NGGLFIMNGS RRRTPGGVFL NLLKNTSPIS EEQIKDIFYI ENQKEYENKK AARKRRTQVL GKKMKQAIKS LNFQEDDDTS RETFASDTNE ALASLDESQE GHAEAKLEAE EAIEVDHSHD

image



Background

Alternative Names	Phosphorylated adapter RNA export protein (RNA U small nuclear RNA
--------------------------	--

Product Name: Recombinant Human PHAX Protein
Catalog #: PEH9424



export adapter protein)

Note

For Research Use Only , Not for Diagnostic Use