


Product Name: Recombinant Human BRD9 Protein
Catalog #: PEH9426



Summary

Name	BRD9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H8M2
Host	E.coli
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DKPLEKPLKL VLKVGGEVST ELSGSGHDSS YYDRSDHER ERHKEKKEKKEK KKKSEKEKHL DDEERRKRKE EKRRREREH CDTEGEADDF DPGKKVEVEP PPDRPVRRACR TQPAENESTP IQQLEHFLR QLQRKDPHGF FAFPVTDAIA PGYSMIKHP MDFGTMKDKI VANEYKSVTE FKADFKLMCD NAMTYNRPDT VYYKLAKKIL HAGFKMMSKQ AALLGNEDTA VEEPVEVVP VQVETAKKSK KPSREVISCM FEPEGNACSL TDSTAEHV L ALVEHADEA RDRINRFLPG GKMGYLKRNG DGSLLYSVVN TAEPDADEEE THPVDLSSLS SKLLPGFTTL GFKDERRNKV TFLSSATTAL SMQNNSVFGD LKSDMELLY SAYGDETVGQ CALSLQEFVK DAGSYSKKVV DDLDDQITGG DHSRTLFLQK QRRNVPMKPP DEAKVGDTLG DSSSVLEFM SMKSYPDVSV DISMLSSLGK VKKELDPDDS HLNLDETTKL LQDLHEAQAE

image



Background

Product Name: Recombinant Human BRD9 Protein
Catalog #: PEH9426



Alternative Names Bromodomain-containing protein 9 (Rhabdomyosarcoma antigen MU-RMS-40.8)

Note

For Research Use Only , Not for Diagnostic Use