


Product Name: Recombinant Human ACTR5 Protein
Catalog #: PEH9428



Summary

Name	ARP5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H9F9
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDSPEELQSY IQKLSIAVEQ AKQKILQAEV NLEVDVVDSK PETPDLEQLE PSLEDVESMN DFDPLFSEET PGVEKPVTTV QPVFNLAAYH QLFVGTERIR APEIIFQPSL IGEEQAGIAE TLQYILDRYP KDIQEMLVQN VFLTGGNTMY PGMKARMEKE LLEMRFPRSS FQVQLASNPV LDAWYGARDW ALNHLDDNEV WITRKEYEEK

image



Background

Alternative Names Actin-related protein 5 (hARP5) (Sarcoma antigen NY-SAR-16)

Note

Product Name: Recombinant Human ACTR5 Protein
Catalog #: PEH9428



For Research Use Only , Not for Diagnostic Use