


Product Name: Recombinant Human UNKL Protein
Catalog #: PEH9429



Summary

Name	UNKL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H9P5
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RGSEQDSKQN HLAVF A AVHP P A PSVSSVA S S LASSAGSG SSSPTALPAP PARALPLGPA SSTVEAVLGS ALDLHLSNVN IASLEKDL E E QDGHD L GAAG PRSLAGSAPV AIPGSLPRAP SLHSPSSAST SPLGSLSQPL PGPV G SSAMT PPQQPPPLRS EPGTLGSAAS SYSPLGLNGV PGSIWDFVSG SFSPSPSPIL SAGPPSSSSA SPNGAELARV

image



Background

Alternative Names Putative E3 ubiquitin-protein ligase UNKL (EC 2.3.2.-) (RING finger protein unkempt-like) (Zinc finger CCCH domain-containing protein 5-like)

Note

Product Name: Recombinant Human UNKL Protein
Catalog #: PEH9429



For Research Use Only , Not for Diagnostic Use