

Product Name: Recombinant Human PLEKHA5 Protein
Catalog #: PEH9431



Summary

Name	PKHA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HAU0
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERPRSAVEQL CLAESTRPRM TVEEQMERIR RHQQACLREK KKGLNVIGAS DQSPLQSPSN LRDNPFRTTQ TRRRDDKELD TAIRENDVKP DHETPATEIV QLKETEPQNV DFSKELKTE NISYEMLFEP EPNGVNSVEM MDKERNKDKM PEDVTFSPQD ETQTANHKPE EHPEENTKNS VDEQEETVIS YESTPEVSRG NQTMAVKSLs

image



Background

Alternative Names	Pleckstrin homology domain-containing family A member 5 (PH domain-containing family A member 5) (Phosphoinositol 3-phosphate-binding protein 2) (PEPP-2)
--------------------------	---

Product Name: Recombinant Human PLEKHA5 Protein
Catalog #: PEH9431



Note

For Research Use Only , Not for Diagnostic Use