


Product Name: Recombinant Human GPN1 Protein
Catalog #: PEH9432



Summary

Name	GPN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HCN4
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LYACSIYKT KLPFIVVMNK TDIIDHSFAV EWMQDFEAFQ DALNQETTYV SNLTRMSLV LDEFYSSLRV VGSAVLGTG LDEFVQVTS AAEEYEREYR PEYERLKKSL ANAESQQQRE QLERLRKDMG SVALDAGTAK DSLSPVLHPS DLILTRGTLT EEDEEADSDT DDIDHRVTEE SHEEPAFQNF

image



Background

Alternative Names	GPN-loop GTPase 1 (EC 3.6.5.-) (MBD2-interacting protein) (MBDin) (RNAPII-associated protein 4) (XPA-binding protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GPN1 Protein
Catalog #: PEH9432

