


Product Name: Recombinant Human XPNPEP3 protein
Catalog #: PEH9442



Summary

Name	XPP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQH7
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIERMQIAGK LTSQAFIETM FTSKAPVEEA FLYAKFEFEC RARGADILAY PPVVAGGNRS NTLHYVKNNQ LIKDGEMLVLL DGGCESSCYV SDITRTWPVN GRFTAPQAEI YEAVLEIQRD CLALCFPGTS LENIYSMMLT LIGQKLKDLG IMKNIKENNA FKAARKYCPH HVGHYLGMDV HDTPDMPRSL PLQPGMVITI EPGIYIPEDD KDAPEKFRGL GVRIEDDVVV TQDSPLILSA

image



Background

Alternative Names	Xaa-Pro aminopeptidase 3 (X-Pro aminopeptidase 3) (EC 3.4.11.9) (Aminopeptidase P3) (APP3)
--------------------------	--

Note

Product Name: Recombinant Human XPNPEP3 protein
Catalog #: PEH9442



For Research Use Only , Not for Diagnostic Use