


Product Name: Recombinant Human BIN3 Protein
Catalog #: PEH9444



Summary

Name	BIN3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQY0
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PKKQIVPKTV ERDFEREYGK LQQLLEEQTTR LQKDMKKSTD ADLAMSKSAV KISLDLLSNP LCEQDQDLLN MVTALDTAMK RMDAFNQEKV NQIQKTVIEP LKKFGSVFPS LNMAVKRREQ ALQDYRRLLQA KVEKYEEKEK TGPVLAKLHQ AREELRPVRE DFEAKNRQLL EEMPRFYGSR LDYFQPSFES LIRAQVVYYS EMHKIFGDLS HQLDQPGHSD EQRERENEAK

image



Background

Alternative Names Bridging integrator 3

Note

Product Name: Recombinant Human BIN3 Protein
Catalog #: PEH9444



For Research Use Only , Not for Diagnostic Use