


Product Name: Recombinant Human STARD7 Protein
Catalog #: PEH9445



Summary

Name	STAR7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQZ5
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QEEELQRSIN EMKRLEEMSN MFQSSGVQHH PPEPKAQTEG NEDSEGKEQR WEMVMDKKHF KLWRRPITGT HLYQYRVFGT YTDVTPRQFF NVQLDTEYRK KWDALVIKLE VIERDVVSGS EVLHWVTHFP YPMYSRDYVY VRRYSVDQEN NMMVLVSRAV EHPSVPESPE FVRVRSYESQ MVIRPHKSFD ENGFDYLLTY SDNPQTVFPR YCVSWMVSSG MPDFLEKLHM ATLKAKNMEI KVKDYISAKP LEMSSEAKAT SQSSERKNEG SCGPARIEYA

image



Background

Alternative Names	StAR-related lipid transfer protein 7, mitochondrial (Gestational trophoblastic tumor protein 1) (START domain-containing protein 7) (StARD7)
--------------------------	---

Product Name: Recombinant Human STARD7 Protein
Catalog #: PEH9445



Note

For Research Use Only , Not for Diagnostic Use