


Product Name: Recombinant Human TXLNG Protein
Catalog #: PEH9452



Summary

Name	TXLNG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NUQ3
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRGGGAEET EAGRGGRRRS PRQKFEIGTM EEAGICGLGV KADMLCNSQS NDILQHQQSN CGGTSNKHSL EEDEGSDFIT ENRNLVSPAY CTQESREEIP GGEARTDPPD GQQDSECNRN KEKTLGKEVL LLMQALNTLS TPEEKLAALC KKYADLLEES RSVQKQMKIL QKKQAQIVKE KVHLQSEHSK AILARSKLES LCRELQRHNK TLKEENMQQA REEEERRKEA TAHFQITLNE IQAQLEQHDI HNAKLRQENI ELGEKLKKLI EQYALREEHI DKVFKHKELQ QQLVDAKLQQ TTQLIKEADE KHQREREFLI KEATESRHKY EQMKQQEVQL KQQLSLYMDK FEFQTTMAK SNEFFTFRQ EMEKMTKKIK KLEKETIWR TKWENNNKAL LQMAEEKTVR

image



Background

Product Name: Recombinant Human TXLNG Protein
Catalog #: PEH9452



Alternative Names

Gamma-taxilin (Environmental lipopolysaccharide-responding gene protein)
(Factor inhibiting ATF4-mediated transcription) (FIAT) (Lipopolysaccharide-specific response protein 5)

Note

For Research Use Only , Not for Diagnostic Use