


Product Name: Recombinant Human GIMAP4 Protein
Catalog #: PEH9454



Summary

Name	GIMA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NUV9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QELHRVELER EKARIREEYE EKIRKLEDKV EQEKRRKKQME KKLAEQEAHY AVRQQRARTE VESKDGILEL

image



Background

Alternative Names GTPase IMAP family member 4 (Immunity-associated nucleotide 1 protein) (IAN-1) (hIAN1) (Immunity-associated protein 4)

Note

For Research Use Only , Not for Diagnostic Use