


Product Name: Recombinant Human THG1L Protein
Catalog #: PEH9457



Summary

Name	THG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NWX6
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FKRKTNWFKR RASKFMTHVA SQFASSYVFY WRDYFEDQPL LYPPGFDGRV VVYPSNQTLK DYLSWRQADC HINNLYNTVF WALIQQSGLT PVQAQGRLLQG TLAADKNEIL FSEFNINYN N ELPMYRKGT V

image



Background

Alternative Names	Probable tRNA(His) guanylyltransferase (EC 2.7.7.79) (Induced in high glucose-1) (IHG-1) (Interphase cytoplasmic foci protein 45) (tRNA-histidine guanylyltransferase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human THG1L Protein
Catalog #: PEH9457

