


Product Name: Recombinant Human THG1L Protein
Catalog #: PEH9458



Summary

Name	THG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9NWX6
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKSKFEYVRD FEADDTCLAH CWVVVRLDGR NFHRFAEKHN FAKPNSRAL QLMTKCAQTV MEELEDIVIA YGQSDEYSFV FKRKTNWFKR RASKFMTHVA SQFASSYVFY WRDYFEDQPL LYPPGFDGRV VVYPSNQTLK DYLSWRQADC HINNLYNTVF WALIQQSGLT PVQAQGRLQG TLAADKNEIL FSEFNINYN ELPMYRKGTV LIWQKVDEVM TKEIKLPTEM EGKKMAVTRT

image



Background

Alternative Names	Probable tRNA(His) guanylyltransferase (EC 2.7.7.79) (Induced in high glucose-1) (IHG-1) (Interphase cytoplasmic foci protein 45) (tRNA-histidine guanylyltransferase)
--------------------------	--

Product Name: Recombinant Human THG1L Protein
Catalog #: PEH9458



Note

For Research Use Only , Not for Diagnostic Use