


Product Name: Recombinant Human KLHL14 Protein
Catalog #: PEH9466



Summary

Name	KLH14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2G3
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPSQQQPPP QEEPSTPSS PDKLLTSPR AINNVLQGC SSIGLRLVLE YLYTANVTLS LDTVEEVLSV SKILHIPQVT KLCVQFLNDQ ISVQNYKQVC KIAALHGLEE TKKLANKYLV EDVLLLNFEED MRALLDSLPP

image



Background

Alternative Names	Kelch-like protein 14 (Protein interactor of Torsin-1A) (Printor) (Protein interactor of torsinA)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use