


Product Name: Recombinant Human STARD9 Protein
Catalog #: PEH9468



Summary

Name	STAR9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2P6
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	INVSLEKDMP GESAVSLKSR SVDRRVSSPV MVAQGGGPTP KWEGKNETGL LEKGLRPKDS SEEFKLPGTK PAYERFQLVA CPQERNPSEC KSQEMLNPNR EPGKKQNKR VNNTDEMARL IRSVMQLENG ILEIESKQNK QVHASHTPGT DKELVFQDQK EQEKTDHAFR PDSSGNPLPS KDQPSSPRQT DDTVFRDSEA GAMEVNSIGN HPQVQKITPN PFRSREGVRE SEPVREHHP AGSDRPARDI CDSLKHTTC REFTNTSLHP QRMKALARAL PLQPRLERSS KNNGQFVKAS ASLKGQPWGL GSLEELETVK GFQESQVAEH VSSSNQEEPQ AQGKVEEMPM QRGGSLSQEEEN

image



Background

Alternative Names	StAR-related lipid transfer protein 9 (START domain-containing protein 9)
--------------------------	---

Product Name: Recombinant Human STARD9 Protein
Catalog #: PEH9468



(StARD9)

Note

For Research Use Only , Not for Diagnostic Use