


Summary

Name	CPNE7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBL6
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTRPLLLKFG RNAGKSTITV IAEDISGNG YVELSFRARK LDDKDLFSKS DPFLELYRVN DDQGLQLVYR TEVVKNNLNP VWEAFKVLSL SLCSCEETRP LKCLVWDYDS RGKHDFIGEF STTFEEMQKA FEEGQAQWDC VNPKYKQKRR SYKNSGVVVL ADLKFHRVYS FLDYIMGGCQ IHFTVAIDFT ASNGDPRNSC SLHYINPYQP NEYLKALVSV GEICQDYDSD KRFSALGFGA RIPPKYEVSH DFAINFNPED DECEGIQGVV EAYQNCLPRV QLYGPTNVAP IISKVARVAA AEESTGKASQ YYILLILTDG VVTDMADTRE AIVRASRLPM SIIIVGVGNA DFTDMQVLDG DDGVLRSRPG EPALRDIVQF VPFRELKNAS PAALAKCVLA

image



Background

Alternative Names Copine-7 (Copine VII)

Product Name: Recombinant Human CPNE7 Protein
Catalog #: PEH9470



Note

For Research Use Only , Not for Diagnostic Use