



Summary

Name	DJB11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBS4
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GVPRASIKD IKKAYRKLAL QLHPDRNPDD PQAQEKFDL GAAYEVLSDS EK RKQYDTYG EEGLKDGHS SHGDIFSHFF GDFGFMFGGT PRQQDRNIPR GSDIIVDLEV TLEEVYAGNF VEVVRNKPVA RQAPGKRRCN CRQEMRTTQL GPGRFQMTQE VVCDECPNVK LVNEERTLEV EIEPGVRDGM EYPPFIGEGEP HVDGEPGDLR FRIKVVKHPI FERRGDDLYT NVTISLVESL VGFEMDITHL DGHKVHISR D KITRPGAKLW KKGEGLPNFD NNNIKGSLII TFDVDFPKEQ LTEEAREGIK QLLKQGSVQK

image



SDS



WB

Background

Alternative Names	DnaJ homolog subfamily B member 11 (APOBEC1-binding protein 2) (ABBP-2) (DnaJ protein homolog 9) (ER-associated DNAJ) (ER-associated
--------------------------	--

Product Name: Recombinant Human DNAJB11 Protein
Catalog #: PEH9472



Hsp40 co-chaperone) (Endoplasmic reticulum DNA J domain-containing protein 3) (ER-resident protein ERdj3) (ERdj3) (ERj3p) (HEDJ) (Human DnaJ protein 9) (hDj-9) (PWP1-interacting protein 4)

Note

For Research Use Only , Not for Diagnostic Use