


Product Name: Recombinant Human CCDC39 Protein
Catalog #: PEH9474



Summary

Name	CCD39_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UFE4
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NEENKLEDQ LSKLKDERAS LQDELREYEE RINSMTSHFK NVKQELSITQ SLCKARERET ESEEHFKAIA QRELGRVKDE IQRLENEMAS ILEKKSDDNK GIFKATQKLD GLKCQMNWDQ QALEAWLEES AHKSDALTL QKYAQQDDNK IRALTLQLER LTLECNQKRK ILDNELTETI SAQLELDKAA QDFRKIHNER QELIKQWENT IEQMQRDGD IDNCALELAR IKQETREKEN LVKEKIKFLE SEIGNNTEFE KRISVADRKL LKCRATYQDH ETSRIQLKGE LDSLKATVNR TSSDLEALRK

image



Background

Alternative Names Coiled-coil domain-containing protein 39

Product Name: Recombinant Human CCDC39 Protein
Catalog #: PEH9474



Note

For Research Use Only , Not for Diagnostic Use