


Product Name: Recombinant Human PACSIN3 Protein
Catalog #: PEH9477



Summary

Name	PACN3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKS6
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EWSLDTQRTI SRKEKGGRSP DEVTLSIVP TRDGTAPPPQ SPGSPGTGQD EEWSDEESPR KAATGVRVRA LYDYAGQEAD ELSFRAGEEL LKMSEDEEQG WCQGQLQSGR IGLYPANYVE

image



Background

Alternative Names	Protein kinase C and casein kinase substrate in neurons protein 3 (SH3 domain-containing protein 6511)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use