


Product Name: Recombinant Human PALD1 Protein
Catalog #: PEH9480



Summary

Name	PALD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULE6
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKAHLSEPPP GKEGPLTYRF QTCLTMQEVF SQHRRACPGL TYHRIPMPDF CAPREEDFDQ LLEALRAALS KDPGTGFVFS CLSGQGRTTT AMVVAVLAFW HIQGFPEVGE EELVSVPDAK FTKGEFQVVM KVVQLLPDGH RVKKEVDAAL DTVSETMTPM HYHLREIIC TYRQAKAAKE AQEMRRLQLR SLQYLERYVC LILFNAYLHL EKADSWQRPF STWMQEVASK AGIYEILNEL GFPELESGED

image



Background

Alternative Names Paladin

Note

Product Name: Recombinant Human PALD1 Protein
Catalog #: PEH9480



For Research Use Only , Not for Diagnostic Use