


Product Name: Recombinant Human CBLC Protein
Catalog #: PEH9481



Summary

Name	CBLC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULV8
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VKIEPCGHLL CSCCLAAWQH SDSQTCPCFR CEIKGWEAVS IYQFHGQATA EDSGNSSDQE GRELELGQVP LSAPPLPPRP DLPPRKPRNA QPKVRLLKGN SPPAALGPQD

image



Background

Alternative Names E3 ubiquitin-protein ligase CBL-C (EC 2.3.2.27) (RING finger protein 57) (RING-type E3 ubiquitin transferase CBL-C) (SH3-binding protein CBL-3) (SH3-binding protein CBL-C) (Signal transduction protein CBL-C)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CBLC Protein
Catalog #: PEH9481

