


Product Name: Recombinant Human ARHGAP26 Protein
Catalog #: PEH9484



Summary

Name	RHG26_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UNA1
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LHLSRKKSSD SKPPSCSERP LTLFHTVQST EKQEQRNSII NSSLESVSSN PNSILNSSSS LQPNMNSSDP DLAVVKPTRP NSLPPNPSPT SPLSPSWPMF SAPSSPMPTS STSSDSSPVR

image



Background

Alternative Names	Rho GTPase-activating protein 26 (GTPase regulator associated with focal adhesion kinase) (Oligophrenin-1-like protein) (Rho-type GTPase-activating protein 26)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ARHGAP26 Protein
Catalog #: PEH9484

