


Product Name: Recombinant Human NDUFAF1 Protein
Catalog #: PEH9495



Summary

Name	CIA30_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y375
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQRKTEGDLQ GDHQKEVALD ITSSEKPDV SFDKAIRDEA IYHFRLKDE IVDHWRGPEG HPLHEVLLEQ AKVWVWQFRGK EDLDKWTVT DKTIGGRSEV FLKMGKNNQS ALLYGTLSE

image



Background

Alternative Names	Complex I intermediate-associated protein 30, mitochondrial (NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use