


Product Name: Recombinant Human GRIP1 Protein
Catalog #: PEH9498



Summary

Name	GRIP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y3R0
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LAGQVVHTET TEVVLTA DPV TGFGIQLQGS VFATETLSSP PLISYIEADS PAERCGVLQI GDRVMAINGI PTEDSTFEEA SQLLRDSSIT SKVTLEIEFD VAESVIPSSG TFHVKLPKKH NVELGITISS PSSRKPGDPL VISDIKKGSV AHRTGTLELG DKLLAIDNIR LDNCSMEDAV QILQQCEDLV KLKIRKDEDN SDEQESSGAI

image



Background

Alternative Names Glutamate receptor-interacting protein 1 (GRIP-1)

Note

Product Name: Recombinant Human GRIP1 Protein
Catalog #: PEH9498



For Research Use Only , Not for Diagnostic Use