


Product Name: Recombinant Human XRN1 Protein
Catalog #: PEH9500



Summary

Name	XRN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZH2
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGGYINESGH LNLPRFEKYL VKLSDFDREH FSEVFVDLKW FESKVGNKYL NEAAGVAAEE ARNYKEKKKL KGQENSLCWT ALDKNEGEMI TSKDNLEDET EDDDLFETEF RQYKRYYMT KMGVDVVSDD FLADQAACYV

image



Background

Alternative Names 5'-3' exoribonuclease 1 (EC 3.1.13.-) (Strand-exchange protein 1 homolog)

Note

For Research Use Only , Not for Diagnostic Use