


Product Name: Recombinant Human ANP32E Protein
Catalog #: PEH9501



Summary

Name	AN32E_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BTT0
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRNRSPEEVT ELVLDNCLCV NGEIEGLNDT FKELEFLSMA NVELSSLARL PSLNKLRKLE LSDNIISGGL EVLAEKCPNL TYLNLSGNKI KDLSTVEALQ NLKLNKSLDL FNCEITNLED YRESIFELLQ QITYLDGFDQ EDNEAPDSEE

image



Background

Alternative Names Acidic leucine-rich nuclear phosphoprotein 32 family member E (LANP-like protein) (LANP-L)

Note

For Research Use Only , Not for Diagnostic Use