


Product Name: Recombinant Human ACTR10 Protein
Catalog #: PEH9502



Summary

Name	ARP10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NZ32
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APSHLMALLT LGINSAMVLD CGYRESLVLP IYEGIPVLNC W GALPLGGKA LHKELETQLL EQCTVDTSVA KEQSLPSVMG SVPEGVLEDI KARTCFVSDL KRGLKIQA AK FNIDGNNERP SPPPNVDYPL DGEKILHILG SIRDSVVEIL FEQDNEEQSV ATLILDSLIQ CPIDTRKQLA ENLVVIGGTS MLPGFLHRL AEIRYLVEKP KYKKALGTKT FRIHTPPAKA NCVAWLGGAI FGALQDILGS RSVSKEYYNQ TGRIPDWCSL NNPPLEMMFD VGKTQPPLMK

image



Background

Alternative Names Actin-related protein 10 (Actin-related protein 11) (hARP11)

Product Name: Recombinant Human ACTR10 Protein
Catalog #: PEH9502



Note

For Research Use Only , Not for Diagnostic Use