



**Product Name: Recombinant Human AK5 Protein**  
**Catalog #: PEH9505**



## Summary

<b>Name</b>	KAD5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9Y6K8
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	35 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	ELAPQETTIT EIKQKLMQIP DEEGIVIDGF PRDVAQALSF EDQICTPDLV VFLACANQRL KERLLKRAEQ QGRPDDNVKA TQRRLMNFKQ NAAPLVKYFQ EKGLIMTFDA DRDEDEVFYD ISMAVDNKLK PNKEAAAGSS DLDPSMILDT GEIIDTGSDY EDQGDDQLNV FGEDTMGGFM EDLRKCKIIF IIGGPGSGKG TQCEKLVEKY GFTHLSTGEL LREELASESE RSKLIRDIME RGDLVPSGIV LELLKEAMVA SLGDTRGFLI DGYPREVKQG EEFGRRIGDP

## image



## Background

**Alternative Names** Adenylate kinase isoenzyme 5 (AK 5) (EC 2.7.4.3) (EC 2.7.4.6) (ATP-AMP transphosphorylase 5)

**Product Name: Recombinant Human AK5 Protein**  
**Catalog #: PEH9505**



---

**Note**

For Research Use Only , Not for Diagnostic Use