


Product Name: Recombinant Human AGXT2 Protein
Catalog #: PEH9510



Summary

Name	AGT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BYV1
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKLAALLPEP LKVIFLVNSG SEANELAMLM ARAHSNNIDI ISFRGAYHGC SPYTLGLTNV GTYKMELPGG TGCQPTMCPD VFRGPWGGSH CRDSPVQTIR KCSCAPDCCQ AKDQYIEQFK DTLSTSVAKS IAGFFAEPIQ GVNGVVQYPK GFLKEAFELV RARGGVCIAD EVQTGFGR LG SHFWGFQTHD VLPDIVTMAK GIGNGFMAA VITTPEIAKS LAKCLQHFNT FGGNPMACAI GSAVLEVIKE ENLQENSQEV GTYMMLLFAK LRDEFEIVGD VRGKGLMIGI EMVQDKISCR PLPREEVNQI HEDCKHMGLL VGRGSIFSQT

image



Background

Alternative Names Alanine--glyoxylate aminotransferase 2, mitochondrial (AGT 2) (EC 2.6.1.44)
((R)-3-amino-2-methylpropionate--pyruvate transaminase) (EC 2.6.1.40)

Product Name: Recombinant Human AGXT2 Protein
Catalog #: PEH9510



(Beta-ALAAT II) (Beta-alanine-pyruvate aminotransferase) (D-AIBAT)

Note

For Research Use Only , Not for Diagnostic Use