


**Product Name: Recombinant Human MAN2C1 Protein**  
**Catalog #: PEH9515**



## Summary

<b>Name</b>	MA2C1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-MBP & C-His
<b>Accession #</b>	Q9NTJ4
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	84 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LVHFPPGDSY GMQGSVEEVL KTVANNRDKG RANHSAFLFG FGDGGGGPTQ TMLDRLKRLS NTDGLPRVQL SSPRQLFSAL ESDSEQLCTW VGELFLELHN GTYTTHAQIK KGNRECEERIL HDVELLSSLA LARSAQFLYP AAQLQHLWRL LLNQFHHDVV TGSCIQMVAE EAMCHYEDIR SHGNTLLSAA AAALCAGEPG PEGLLIVNTL PWKRIEVMAL PKPGGAHSLA LVTVPSMGYA PVPPPTSLQP LLPQQPVFVV QETDGSVTLD NGIIRVKLDP TGRLTSLVLV ASGREAIAEG AVGNQFVLFV DVPLYWDAWD VMDYHLETRK PVLGQAGTLA VGTEGGLRGS AWFLLQISPN SRLSQEVVLD

## image



## Background

**Alternative Names** Alpha-mannosidase 2C1 (EC 3.2.1.24) (Alpha mannosidase 6A8B) (Alpha-D-

**Product Name: Recombinant Human MAN2C1 Protein**  
**Catalog #: PEH9515**



---

mannoside mannohydrolase) (Mannosidase alpha class 2C member 1)

**Note**

For Research Use Only , Not for Diagnostic Use