


Product Name: Recombinant Human ANKRD34A Protein
Catalog #: PEH9521



Summary

Name	AN34A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q69YU3
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GASALVHALD RGDRETLATL LDACKAKGTE VIIIITDTSP SGTKKTRQYL NSPPSPGVED PAPASPSPGF CTSPSEIQLQ TAGGGGRGML SPRAQEEEEK RDVFEFPLPK PPDDPSPSEP LPKPPRHPPK PLKRLNSEPW GLVAPPQPVP PTEGRPGIER LTAEFNGLTL TGRPRLSRRH STEGPEDPPP WAEKVTSGGP LSRRNTAPEA QESGPPSGLR QKLSRMEPVE LDTPGHLCPD SPESSRLSLE RRRYSASPLT LPPAGSAPSP

image



Background

Alternative Names Ankyrin repeat domain-containing protein 34A

Product Name: Recombinant Human ANKRD34A Protein
Catalog #: PEH9521



Note

For Research Use Only , Not for Diagnostic Use